
















BASIC ELECTRONIC AIDS TO DAILY LIVING – Battery and Simple Electrical Devices

NAME	PHOTO	CONTROL MODE			SWITCH INPUT	SWITCH OUTPUT	COST	PROD. #	DIMENSIONS	COMMENTS	
		Direct	Latch	Time sec							Time min
Battery Transmission											
Single SLAT Ablenet			X	0 - 60	0 - 60	1	1	\$75	100-SSLAT	2.25 x 3.75 x 1	
Choice SLAT Ablenet			X	0 - 60	0 - 60	2	2	\$95	100-CSLAT	2.25 x 3.75 x 1	2nd device will not activate until 1st device stops
Dual SLAT Ablenet			X	0 - 60	0 - 60	2	2	\$95	100-DSLAT	2.25 x 3.75 x 1	allows 2 people to use at one time
Tash Switch Latch Ablenet		X	X	2 - 52		2	2	\$125	43880	?	
LinkSwitch Adaptivation		X	X	1-60	1-60	1-2	1-2	\$106	DLT-202	3 x 4.5 x 2.5	alarm mode
Switch Scanner Adaptivation		X	X	1-60	1-60	1-2	1-6	\$159	SWSC	5.75x3.75 x 1.25	scans up to 6 outputs
5-in-1 Switch Modifier Enabling Devices		X	X	0-120		1	1	\$89.95	566	6.75 x 2.5 x 1.5	intention mode (adjusts activation time). Includes adj. pressure switch
Switch Modifier Enabling Devices			X	0 - 120		1	1	\$59.95	605	4.5 x 2.5 x 1	

Ablenet, inc. 800-322-0956 www.ablenetinc.com
 Adaptivation 800-723-2783 www.adaptivation.com
 Enabling Devices/Toys for Special Children 800-832-8697 www.enablingdevices.com

BASIC ELECTRONIC AIDS TO DAILY LIVING – Battery and Simple Electrical Devices

NAME	PHOTO	CONTROL MODE				SWITCH INPUT	SWITCH OUTPUT	COST	PROD. #	DIMENSIONS	COMMENTS
		Direct	Latch	Time sec	Time min						
Electrical Transmission											
PowerLink 3 Control Unit Ablenet		X	X	0-60	0-60	2	4	\$189	100-PL3D	9 x 5 ½ x 2 ¼	up to 1700 watts total remote switch access with AirLink (\$89)
Electralink Ablenet		X	X	1-60	1-60	1	1	\$200	34410	?	up to 1500 watts, 15A. Can be used with wireless Ultra One Transmitter
FreeHand Adaptivation		X	X	1-60	1-60	1	1	\$115	FH-ME	3.5 x 3.5 x 3.5	with transceiver module \$145
FreeSwitch Max Adaptivation		X	X	1-60	1-60	1-4	1-8	\$299	FSM-100	4 x 6.5 x 1.8	visual and auditory scanning (speech) X10 modules required Direct access option
Wireless Pal Pad Adaptivation		X	X			1	1	\$95	WPP	3.5 x 4	with transceiver module \$125
1-Switch Single- Appliance Unit Enabling Devices		X	X			1	1	\$107.95	1490	3 x 2.5 x 2	remote up to 600 watts
Environmental Controls Enabling Devices			X			1	4	\$149.95	591	4.5 x 2.5 x 1	Includes 2 receivers. 2 more \$63.95.

Basic EADLs provide direct, latched or timed on/off switch control of battery operated toys or simple electrical appliances. For control of more devices and features, please refer to the Multifunction EADLs chart.

Ablenet, inc.	800-322-0956	www.ablenetinc.com
Adaptivation	800-723-2783	www.adaptivation.com
Enabling Devices/Toys for Special Children	800-832-8697	www.enablingdevices.com