







**BASIC ELECTRONIC AIDS TO DAILY LIVING – Battery and Simple Electrical Devices**

NAME	PHOTO	CONTROL MODE				SWITCH INPUT	SWITCH OUTPUT	COST	PROD. #	DIMENSIONS	COMMENTS
		Direct	Latch	Time sec	Time min						
<b>Battery Transmission</b>											
Single SLAT AbleNet			X	0 - 60	0 - 60	1	1	\$99	100-SSLAT	2.25 x 3.75 x 1	
Choice SLAT AbleNet			X	0 - 60	0 - 60	2	2	\$119	100-CSLAT	2.25 x 3.75 x 1	2nd device will not activate until 1st device stops
Dual SLAT AbleNet			X	0 - 60	0 - 60	2	2	\$119	100-DSLAT	2.25 x 3.75 x 1	allows 2 people to use at one time
Jelly Beamer with SLAT Big Beamer with SLAT AbleNet		X	X	0-60	0-60	1	1	\$208	100-34000 100-33700		Wireless switch. Has external switch jack, as well.
LinkSwitch Adaptation		X	X	1-60	1-60	1-2	1-2	\$1015	DLT-202	3 x 4.5 x 2.5	alarm mode
Switch Scanner Adaptation		X	X	1-60	1-60	1-2	1-6	\$159	SWSC	5.75x3.75 x 1.25	scans up to 6 outputs
5-in-1 Switch Modifier Enabling Devices		X	X	0-120		1	1	\$64.95	566	6.75 x 2.5 x 1.5	intention mode (adjusts activation time). Includes adj. pressure switch
Switch Modifier Enabling Devices			X	0 - 120		1	1	\$64.95	605	4.5 x 2.5 x 1	

AbleNet, inc. 800-322-0956 www.ablenetinc.com  
 Adaptation 800-723-2783 www.adaptation.com  
 Enabling Devices/Toys for Special Children 800-832-8697 www.enablingdevices.com

**BASIC ELECTRONIC AIDS TO DAILY LIVING – Battery and Simple Electrical Devices**

NAME	PHOTO	CONTROL MODE		Time sec	Time min	SWITCH INPUT	SWITCH OUTPUT	COST	PROD. #	DIMENSIONS	COMMENTS
		Direct	Latch								

<b>Electrical Transmission</b>											
PowerLink 4 Control Unit AbleNet		X	X	0-60	0-60	2	4	\$229	100-10701	9 x 5 1/2 x 2 1/4	up to 1700 watts total remote switch access with Jelly Beamer
FreeHand Adaptation		X	X	1-60	1-60	1	1	\$115	FH-ME	3.5 x 3.5 x 3.5	with transceiver module \$145
FreeSwitch Max Adaptation		X	X	1-60	1-60	1-4	1-8	\$299	FSM-100	4 x 6.5 x 1.8	visual and auditory scanning (speech) X10 modules required Direct access option
Wireless Pal Pad Adaptation		X	X			1	1	\$95	WPP	3.5 x 4	with transceiver module \$125
Continuous Control Unit Enabling Devices		X	X			1	1	\$112.95	1490	3 x 2.5 x 2	remote up to 600 watts
Environmental Controls Enabling Devices			X			1	4	\$152.95	591	4.5 x 2.5 x 1	Includes 2 receivers. 2 more \$67.95.

Basic EADLs provide direct, latched or timed on/off switch control of battery operated toys or simple electrical appliances. For control of more devices and features, please refer to the Multifunction EADLs chart.

AbleNet, inc. 800-322-0956 www.ablenetinc.com  
 Adaptation 800-723-2783 www.adaptation.com  
 Enabling Devices/Toys for Special Children 800-832-8697 www.enablingdevices.com