




**BASIC ELECTRONIC AIDS TO DAILY LIVING – Infrared Control**

NAME	PHOTO	SWITCH INPUT	SWITCH OUTPUT	COST	PROD#	DIMENSIONS	COMMENTS
<b>Infrared Transmission</b>							
a. Mini Relax - Scanning b. Mini Relax with Jacks c. Mini Relax with X-10 Ablenet		a. 1 b. 1-6 c. 1	a.6 IR cmds b.6 IR cmds, c.6 IR cmds 1 X-10	\$300 \$325 \$400	82050 82060 82070	6 x 3.5 x 1.5	a. auditory and visual feedback, auto scan b. one switch required for each function c. includes X-10 radio transceiver
Wireless TV Remote Enabling Devices		direct or 5 switches	TV	\$142.95	1521	15.75 x 10.5 x 3.5	2" buttons or any 5 switches to control power, channels, volume
TV Remote Module Enabling Devices		1-5 switches	TV	\$72.95	5150	10 x 4 x 2.5	switch for each desired function: power, volume up, volume down, channel up, channel down

Basic EADLs send a limited amount of infrared signals to a limited amount of devices. For control of more devices and features, please refer to the Multifunction EADLs chart.

Ablenet, Inc. 800-322-0956 [www.ablenetinc.com](http://www.ablenetinc.com)  
 Enabling Devices/Toys for Special Children 800-832-8697 [www.enablingdevices.com](http://www.enablingdevices.com)

Michelle L. Lange, OTR, ABDA, ATP, Access to Independence, [MichelleLange@msn.com](mailto:MichelleLange@msn.com), 5/08.