






## Infrared Controlled Telephones

An infrared controlled telephone can be used by virtually any Electronic Daily Living Aid (EADL, formally Environmental Controls) that can send infrared signal. Here is a comparison of features:

Features	<b>GewaTel 200</b> ZYGO Industries, Inc. 800-234-6006 <a href="http://www.zygo-usa.com">www.zygo-usa.com</a>	<b>Infrared Home or Office Accessible Landline Telephone</b> Broadened Horizons 612-851-1040 <a href="http://www.broadenedhorizons.com">www.broadenedhorizons.com</a>	<b>Prism Prism D</b> Ablenet, Inc. 800-322-0956 <a href="http://www.ablenetinc.com">www.ablenetinc.com</a>	<b>SiPhone Quickphone Quickphone f</b> Dr. Hein <a href="http://www.SiCare.de">www.SiCare.de</a> (Germany)	<b>TS Phone with TS Controller</b> Technical Solutions 03 9737 9000 <a href="http://www.tecsol.com.au">www.tecsol.com.au</a> (Australia)
Photo					
cost	\$1,380.00	\$699.00	To be determined	?	\$360 (GST) phone only
access	-IR signal -Keypad	-IR signal -Keypad -Voice with voice dialer	-IR signal -Keypad	-IR signal -Keypad	-Single switch for answer and hang-up only -IR signal with TS Controller -Keypad
user can build #s	Yes	Yes	Yes	Yes	Yes
prestored #s held	50	50	5 speed dial Prism relies on directory in SGD Prism D adds directory of 100	SiPhone 200 Quickphone 50	20
redial	Yes	Yes	Yes	Yes	Yes
call waiting	Yes	Yes	Yes	Yes	Yes
battery back-up	Yes	Yes	Yes	Yes	Yes
Display	Yes	Yes	Yes	Yes	Yes
Volume Control	Yes	Yes	Yes	Yes	Yes
Privacy headset	\$115 earphone \$340 microphone	Wired privacy headset \$29.95 Wireless headset \$349.95 Expansion Microphone \$149.95	Pending	Optional	2 Options \$297 or \$498 (GST)

Interfaces	-any learning IR transmitter -any IR transmitting EADL -GEWA Prog	-any learning IR transmitter -any IR transmitting EADL -VoiceIR Voice Controller	-any learning IR transmitter -any IR transmitting EADL -SiCare Standard, Light II, Relax III	-any learning IR transmitter -any IR transmitting EADL -SiCare Standard	-any learning IR transmitter -any IR transmitting EADL
Transmission	Analog landline	Analog landline	Analog landline	Analog ISDN	?
Other	-can be used as standard phone. -commands pre-programmed into Quickie iQ electronics.	-can be used as a standard phone -voice dialer \$249.95 -same base phone as GewaTel 200	-same base phone as Quickphone -Prism D	SiPhone is cordless SiPhone takes pictures, sends via MMS "f" version has emergency call function, requires remote button press -can be used as a standard phone	-can be used as a standard phone

Michelle L. Lange, OTR, ABDA, ATP, Access to Independence, [MichelleLange@msn.com](mailto:MichelleLange@msn.com), 5/08.